

PIK3C2A-phosphoinositide-3-kinase, class 2, alpha polypeptide Gene

Reference: PTXBC113658

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPK; PI3-K-C2(ALPHA); PI3-K-C2A; PI3K-C2-alpha; PI3K-C2alpha; phosphatidylinositol 4-phosphate 3-kinase C2 domain-containing subunit alpha; C2-containing phosphatidylinositol kinase; phosphatidylinositol-4-phosphate 3-kinase C2 domain-containing subunit a		

For research use only.

GENE DESCRIPTION	phosphoinositide-3-kinase, class 2, alpha polypeptide
GENE ID	5286
ORIGIN SPECIES	Human
ACCESSIONS	BC113658
VECTOR	pDONR223
NCBI SYMBOL	PIK3C2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr