

**AHSA2-AHA1, activator of heat shock 90kDa protein ATPase homolog 2 (yeast) Gene**  
Reference: PTXBC113671

|                |                                                                                                         |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                     | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | activator of 90 kDa heat shock protein ATPase homolog 2; AHA1, activator of heat shock protein ATPase 2 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | AHA1, activator of heat shock 90kDa protein ATPase homolog 2 (yeast)                              |
| <b>GENE ID</b>          | 130872                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC113671                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | AHSA2                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

**ORF (Only for reference, product provided without STOP codon )**

atgattttaccaacgaagctatggcaaccaggaattgactgtcaaaagaaactgagtggaataccctgcaggttcaggcctcatccagtggcactgggttaaggattcccactgtggtcttcacatgatggaactgtttgacacaactgtagagcagctgtatagcattctactgt  
aaaggagttaacaaataagaatcatcatgaaatggagatgtgggaactggccagaagaacactatgccatggttcactgaattttgtgcctactctaggcacaacagaattacaattgaaggagttcctatctctgttaagaagagaacatgaaattctgttggcagaagcagcattttg  
aagaaataaagggttcactgcagctgacccccctaaatggttga