

LOC339535-hypothetical LOC339535 Gene

Reference: PTXBC113688

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC339535		

For research use only.

GENE DESCRIPTION	hypothetical LOC339535
GENE ID	339535
ORIGIN SPECIES	Human
ACCESSIONS	BC113688
VECTOR	pDONR223
NCBI SYMBOL	LOC339535
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttcccaaaacagtccccctgtcatcaggaagcagtaaacactggtcttcttctacgttattcttaacagcagttatatgtacttcttttagacgggggaataagacagccaggtgggaggggggtctcttgaaaaactccatctggcctgtgcactagggtggaacctcaggaa
gttgatgacattttagctgggagcagctggccccctcttctctgtgtggaacctgaaattccactggttggcaggaagcgtctagcagggaactccggtttgcagagcatccctgtttcccccgttttcccttttccaccaataaaacctgttctactcaccctcaaacctctgaagc
ctaaatttttggctgtgggacggacaaagatgccttagctga