



Reference: PTXBC113709

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EAG; EAG1; Kv10.1; TMBTS; ZLS1; h-eag; hEAG1; potassium voltage-gated channel subfamily H member 1; EAG channel 1; ether-a-go-go potassium channel 1; ether-a-go-go, Drosophila, homolog of; potassium channel, voltage gated eag related subfamily H, member 1		

For research use only.

GENE DESCRIPTION	potassium voltage-gated channel, subfamily H (eag-related), member 1
GENE ID	3756
ORIGIN SPECIES	Human
ACCESSIONS	BC113709
VECTOR	pDONR223
NCBI SYMBOL	KCNH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]