



Reference: PTXBC113714

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LRP-10; LRP9; MST087; MSTP087; low-density lipoprotein receptor-related protein 10; LDL receptor related protein 10		

For research use only.

GENE DESCRIPTION	low density lipoprotein receptor-related protein 10
GENE ID	26020
ORIGIN SPECIES	Human
ACCESSIONS	BC113714
VECTOR	pDONR223
NCBI SYMBOL	LRP10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr