

C6orf25-chromosome 6 open reading frame 25 Gene

Reference: PTXBC113719

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G6b-B; NG31; protein G6b; immunoglobulin receptor; chromosome 6 open reading frame 25		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 25
GENE ID	80739
ORIGIN SPECIES	Human
ACCESSIONS	BC113719
VECTOR	pDONR223
NCBI SYMBOL	C6orf25
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgtttctgcagctgtaccgctgctctcagggcccaagggaacctgggcttctctgacggccgctggggaccgggtgaatctctctcggaggagtctctcatccatccgctgggtctgggcaccagcttccggcctgaaggcctgtccaaaggacgc
cgaccgatctgtgggctcttctgagcgggacccccaccgtgctctccctccagccttctgctggcgctacgtccctggactctgtatccggcgctggagctcctcttgagcgcgggggactcgggcactttttctgcaaggccgcccacgaggacgagccgtacagtgtt
cacgtgctggggacaggacattgcaaggccccgggctaccatgggtccgtgtatccccagctcctgacccgctgctggcgctgggttggctcggactgggagcttgggcctgtctgtgtgctgcacaggcgctgccccgcaaccgattcgaccactccctagatt
gctctgtccccccacatagctccactgtgaaaaccgagccccagaggccagtaagggaagagcccaagattccaggggacctggaccaggaaccggttaaggcatgggatgggaagggtatag