

FAM127C-family with sequence similarity 127, member C Gene

Reference: PTXBC113720

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM127C; CXX1c; MAR8B; mammalian retrotransposon derived protein 8B; family with sequence similarity 127 member C		

For research use only.

GENE DESCRIPTION	family with sequence similarity 127, member C
GENE ID	441518
ORIGIN SPECIES	Human
ACCESSIONS	BC113720
VECTOR	pDONR223
NCBI SYMBOL	FAM127C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaggtcgagtcagctgatgaagccctctggtcgccctccggcccgccgctgctgaggaacccgattcccttcccgagacgtttgatggcgataccgacggctcccgagttcatcgtgcagacgagctctacatgttcgtggacgagaacacgttccaac
gacccctgaaggtgacgttctcatcaccgcctcagggcccgccctgcagtggtgatccctacatcaagaaggagagccccctcctcagtgattaccgggcttctggtgagatgaagcgggtctttgatgggaggagacgaggacttctag