

PRDX5-peroxiredoxin 5 Gene

Reference: PTXBC113723

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACR1; AOEB166; B166; HEL-S-55; PLP; PMP20; PRDX6; PRXV; SBB110; prx-V; peroxiredoxin-5, mitochondrial; Alu co-repressor 1; TPx type VI; antioxidant enzyme B166; epididymis secretory protein Li 55; liver tissue 2D-page spot 71B; peroxiredoxin V; peroxisoma		

For research use only.

GENE DESCRIPTION	peroxiredoxin 5
GENE ID	25824
ORIGIN SPECIES	Human
ACCESSIONS	BC113723
VECTOR	pDONR223
NCBI SYMBOL	PRDX5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggactagctggcgtgtgcgcctgagacgtcagcgggctatatactctcgtggggccggcggtcagctcgggcagcggcagcaagcgggtcagtgaggagagtggtggcgtctggcgggggtccgcagttcagcagagccgctgcagccatggcccaatcaaggtgg
gagatgccatcccagcagtgagggtgttgaaggggagccagggaacaagggtgaacctggcagagctgttcaagggaagaagggtgtgtctgtttggagttctctggggccttcacctggatgttccaagacacacctgccaggggtttgtggagcaggtgaggtctgaagccaa
gggagtcaggtgtgtgcctgtctgagtgtaatgatgcctttgtgactggcagtggtggccgagcccacaaggcgggaaggcaaggttcggctcctggctgatccactggggcctttgggaaggagacagactattactagatgattcgctgtgtccatcttgggaatcagctctca
agaggttccatggtgtacaggtgcatagtgaggccctgaatgtggaaccagatggcacaggcctcacctgcagcctggcacccaatatcatctcacagctctga