

CXorf42-chromosome X open reading frame 42 Gene

Reference: PTXBC113730

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	chromosome X open reading frame 42
GENE ID	158801
ORIGIN SPECIES	Human
ACCESSIONS	BC113730
VECTOR	pDONR223
NCBI SYMBOL	CXorf42
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcccggttttcaggctgcagtagctccggtgtccggttctgcagtgccgagagggagcctcgggctccggcgggggaatgtgttcagttgtccaagagcccaagccagcaaaagccgcccgttcccgctggggcgcaactctactggaactctgttctctgactgggga
atggaaatgggctcagccaaggctccgaaaccagccctgccgtccagttgtgttcgcatctcaagagcggctgactgcgcgtggaccatcactttcgtctgctcttgcgcctattgcggcaactactacgctcctacgggacgaggagagggtgaatcgtggcag
agaaaaatagtagtctag