

CD274-CD274 molecule Gene

Reference: PTXBC113734

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD274 molecule; CD274 antigen; B7-H; B7H1; PD-L1; PDCD1L1; PDCD1LG1; PDL1; programmed cell death 1 ligand 1; B7 homolog 1; PDCD1 ligand 1		

For research use only.

GENE DESCRIPTION	CD274 molecule
GENE ID	29126
ORIGIN SPECIES	Human
ACCESSIONS	BC113734
VECTOR	pDONR223
NCBI SYMBOL	CD274
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggatatttgcctgtttatattcatgacctactggcatttgcgaacgatttactgtcacgggtcccaaggacatatgtgttagagtatgtagcaatgacaattgaatgcaattcccagtagaaaaacaattagacctggctgcactaattgtctattgggaatggagataagaac
attattcaatttgcctgagagaggaagacctgaaggttcagcatagtagctacagacagaggccggctgttgaaggaccagctctccctgggaattgctgcacttcagatcacagatgtgaaatgacaggtgaggggtgaccgtgcatgacgctatgggtgccgactacaa
gcgaattactgtgaaagtcaatgcccatatacaaaaaaacaacaaagaatttgggtgtggtatccagctcctgaacatgaactgacatgtcaggctgaggggtaccccaaggccgaagtcatctggacaagcagtgaccatcaagtcctgagtggtgaagaccaccaccacaaftcca
agagagaggagaagcttttcaatgtgaccgacactgagaatcaacacaacaactaatgagatttttactgcacttttagagattagatcctgaggaaaaccatacagctgaattggatcccgagaactacctctggcacatcccaaatgaaaggactcacttggaattctggagacc
atcttattatgccttggtagcactgacattcatctccgtttaagaaaaggagaatgatgatgtgaaaaatgtggcatccaagatacaaaactcaagaagcaagtgatcacatttggaggagacgtaa