

AES-amino-terminal enhancer of split Gene

Reference: PTXBC113735

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AES-2; AES-1; ESP1; GRG; GRG5; Grg-5; TLE5; amino-terminal enhancer of split; gp130-associated protein GAM; groucho-related protein 5		

For research use only.

GENE DESCRIPTION	amino-terminal enhancer of split
GENE ID	166
ORIGIN SPECIES	Human
ACCESSIONS	BC113735
VECTOR	pDONR223
NCBI SYMBOL	AES
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtcacaagaatggctccctcaggaaagcggcattaccgctgctttctgcagaaaaggaaactaaggctcagcaagaaccaccgccagcagagccaaggtcacagagcagtcggtggcagcgtccaggtcgtgccacagcaggccggcggtctcagcgggcagc
cgggtcccttttttgacagatggggaaccgaggccggcggtcgtccggggctcctcgcacctaccgccagcaactcaaatccaccctcggactcctgcgacccatcaaaagacgaatttcagctactgcaagctcagttaccacagcctcaagctcgaatgtgacaagtggccagt
gagaagtcagagatgcagcgtcactatgtgatgtactacgagatgtctcagcgttgaacatcgagatgcacaaacggctgagatcgtcaaaaggctgaacgggatttgtgccaggtcctgcctactctcccaagagcaccagcagcaggtcttggagccattgagaggccaa
gcaggtcacgcgtcccgagctgaactctatccgacagcagctccaagccaccagctgtccagctgcaggccctggcccttgacccactaccgtgggctgcagccgctcgtccggcggtcagcgcaggcaccggcctctctcgtgtccgcgtgggttc
caggccacctctccaaggaagacaagaacgggcacgatggtgacacccaccaggagatgatggcgagaagtcgattag