

C5-complement component 5 Gene

Reference: PTXBC113738

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	prepro-C5; complement C5; C5D; C5a; C5b; CPAMD4; ECLZB; C3 and PZP-like alpha-2-macroglobulin domain-containing protein 4; C5a anaphylatoxin; anaphylatoxin C5a analog; complement component 5		

For research use only.

GENE DESCRIPTION	complement component 5
GENE ID	727
ORIGIN SPECIES	Human
ACCESSIONS	BC113738
VECTOR	pDONR223
NCBI SYMBOL	C5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr