



Reference: PTXBC113741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RNA-binding protein MEX3D; MEX-3D; MEX3; OK/SW-cl.4; RKHD1; RNF193; RING finger and KH domain-containing protein 1; RING finger protein 193; bcl-2 ARE RNA binding protein; ring finger (C3HC4 type) and KH domain containing 1; ring finger and KH domain cont		

GENE DESCRIPTION	mex-3 homolog D (C. elegans)
GENE ID	399664
ORIGIN SPECIES	Human
ACCESSIONS	BC113741
VECTOR	pDONR223
NCBI SYMBOL	MEX3D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]