

TMEM64-transmembrane protein 64 Gene

Reference: PTXBC113828

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 64		
For research use only.			

GENE DESCRIPTION	transmembrane protein 64
GENE ID	169200
ORIGIN SPECIES	Human
ACCESSIONS	BC113828
VECTOR	pDONR223
NCBI SYMBOL	TMEM64
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctcagctattactgatctcattacccaactatctgatggcatcttgggtggactgcttctaccagcttctgaattcttacttgggtaccacctcgccgacaatggaagatgtcattgcagaacagagtgttagtggatatgtttttgtttacagattattataagtaggcctcatgtttatgtatgtcagctcagctcaagtgggaattgaatgcagctattgtagcttgtgaaatgggaactgaaattctcttgggttaaggcaatcaaccaataaccagtggtcttctattctacaacaagaggacccctaacatttctggagggtggaatcaatgtttgtatga