

## CD3G-CD3g molecule, gamma (CD3-TCR complex) Gene

Reference: PTXBC113830

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD3g molecule; CD3g molecule, gamma (CD3-TCR complex); CD3g molecule, epsilon (CD3-TCR complex); CD3g antigen, gamma polypeptide (TiT3 complex); CD3-GAMMA; IMD17; T3G; T-cell surface glycoprotein CD3 gamma chain; T-cell antigen receptor complex, gamma sub		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	CD3g molecule, gamma (CD3-TCR complex)
<b>GENE ID</b>	917
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC113830
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	CD3G
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggaacaggggaaggcctggctgctcctcatctggctatcattcttcaaggctctggccagtcacaaaggaaccactgggtaaggtgatgactatcaagaagatggctcgcttctgacttgtgatgcagaagccaaaataatcacatggttaagatgggaagatgac  
ggcttctaactgaagataaaaaaaatggaatctgggaagtaatccaaggaccctcgaggatgtatcagtgtaaaaggatcacagaacaagtcaaaaccactccaagtgtattacagaatgtgtcagaactgcatgaaactaatgcagccacatctggcttcttctgaaatcgt  
cagcatttctgcttctgcttgggtctactcattgctggacaggatggagttccagctcgagagcttcagacaagcagactctgtgcccaatgaccagctctaccagccctcaaggatcagaagatgaccagtacagccacctcaaggaaccagttgaggaggaattga