



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

OR10H3-olfactory receptor, family 10, subfamily H, member 3 Gene

Reference: PTXBC113844

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	olfactory receptor 10H3; olfactory receptor OR19-24; olfactory receptor family 10 subfamily H member 3		

For research use only.

GENE DESCRIPTION	olfactory receptor, family 10, subfamily H, member 3
GENE ID	26532
ORIGIN SPECIES	Human
ACCESSIONS	BC113844
VECTOR	pDONR223
NCBI SYMBOL	OR10H3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctctgcagaactacagaacatcatctgaatttaattctctctgcctctcagccctcccccagcagctctcgtctgtctgtctctcgtctacacctctgattgtctctgttcaacattgttggaacctcttctatcatggccacagtttgatgaacgcagactccacaccccatgtacctctctgtgtgcccctcccatctctgagattctgttcaactgttgccatccccctcgcatgctggctgatctgtctctcaccacgtttccatcaacatttttggtctgtgcccattcagatgttcttctctcatgttttgcttccactactcctctcttctcatggctatgggtatgatcaactacgtgacacatctgcacccaactgcattacaacatctgtaattgagctcccctgtgtgcccattctgtgcccctagggctggggctgtggctcgtgcacgtggagatgggtgacaaatgatggttgacaaatgatggttgcacacctcaactctctgctcttaattgtgatccacatttctctgtcatgctcttcccctctgaagtggcctgtggagcagaagacatctctgcatcatgctgtagctggtgtgtgacaccccctgatggctgtttgttctcatcatctctctcttcttctcattatgggtggccatctctggttggccatctctgtagggccggcggcacaagactcttcccactctgtatccacactcactgtgtggctatgactcatatgtttgtctctctctatctcaaaacccaagggccctcatctctgtctgtacagtgtatgctctgtatggccacacatactgtctaccctctctccagcccaatttttcagtctaaaggaaacaaaggagctgaagaaatgccataaataaaaacttttcagaagaagctctcccctcttaagctctcaaa

More info: Please contact us for any question: contact@proteogenix.fr