



Reference: PTXBC113859

|                |                                                                                                                                                                                                                                                                 |                     |                |
|----------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CD45-AP; LPAP; protein tyrosine phosphatase receptor type C-associated protein; CD45 associated protein; PTPRC-associated protein; lymphocyte phosphatase-associated phosphoprotein; protein tyrosine phosphatase, receptor type, c polypeptide-associated prot |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | protein tyrosine phosphatase, receptor type, C-associated protein                                 |
| <b>GENE ID</b>          | 5790                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC113859                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PTPRCAP                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

atggctctgcgccctgcacctaggcctcgggafgctgctggccctccagggggccttgggctcggggtggcagcgcggaggacagcgtgggctccagctcgtfcaaccgtgttctcgtgctgctgctcctactcgtggccactgctgtagcactggccctggcggccctcagccgtgac  
tcagggggctactaccacccggcccgctaggtgcccggctgtgtggggccgcagcggggcctgtcctgtggccagccccagggctgctgctgctcagggcccagctgagctgtgggtccacagacaaatgacctgtgagcgcagaggaggatgagcagggacacagactatgaccagctc  
cgcgatgctgctgcccctcagcctgacccctggcggcggcggcagcgcagcagaatgtgtgagagcggcgcagccgacagcagggctcccgctgctggctgctgaggaaagccagagacagctgacacacaggagggcgacctggctcctgctggccctccagaccagcagcgggggcagctgctgagggcc  
ctctcgaatgacctcagcggcctcttctgcacccgcagcgcctggatgacacgcacggcgcagcggcggccggccggccctccagctccagctgtag