

RNF130-ring finger protein 130 Gene

Reference: PTXBC113864

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF130; G1RZFP; GOLIATH; g1-related zinc finger protein; goliath homolog; ring finger protein 130		

For research use only.

GENE DESCRIPTION	ring finger protein 130
GENE ID	55819
ORIGIN SPECIES	Human
ACCESSIONS	BC113864
VECTOR	pDONR223
NCBI SYMBOL	RNF130
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcaccagcactggagatatttctgtcatgatacagaattgagggtaaggatattttagttatctggagaaaacatctctgtacaatgacaatagctgttggaaactgaatgccaccgaagaactcagccgtggctctctagtctctgtgtcaatcctttattgtttgatga
ttatttctcagcatggctcatattctacttcattcagaagatcaggtacacaaatgcacgcagcaggaaccagcgtctctcggagatgcagccaagaagccatcagtaaatgacaaccaggacagtaaaaggggtgacaaggaaactgacccagactttgatcattgtgcagctgca
tagagagctataagcagaatgatgtcgtccgaattctccctgcaagcatgtttccacaaatcctgctggatccctgcttagtgaacattgtacctgtctatgtcgaactaataattgaaggccctgggaattgtgccgaattgcatgtactgataacgtagcattcgataggaaagg
ctcaccagaaccaagctgttaaccgaagatcagccctggcgacctcggcgacacaccccttggccttgagccacttgaactcgggagatcaccctctctcaggtggggagctcactccgagaacaggagaataacattgcagtaacaaaagaatggtttattatgccagt
ttggcctcctcagtgccctcacactctgtacatgatcatcagagccacagctagcttgaatgctaagtagtagaatggttttga