

UPK2-uropalakin 2 Gene

Reference: PTXBC113900

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UP2; UPII; uropalakin-2; uropalakin II; uropalakin 2		
For research use only.			

GENE DESCRIPTION	uropalakin 2
GENE ID	7379
ORIGIN SPECIES	Human
ACCESSIONS	BC113900
VECTOR	pDONR223
NCBI SYMBOL	UPK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaccctgctgccatccggacctgcccttgatctgattctgctggctctgctgccccaggggctgcagactcaacatctcaagcctctctggctgctgctcccccgcactaacggagagcctgctggtgcttccccctgtcacctcacaggaggcaatgccacactgatgg
tccggagagccaatgacagcaaatgtggtgacgtccagctttgtggtgctccctgctggcgaggggaactggtgagtggtggtgacagtggtgctggttcacagtcactcggctcagtgatcacagggtgacaaacctgtgcccaggaacaaattctacatttctacctaagttaa
gaaggggacagccactgagtcagcagagagatccaatgtccacactcctcgaaggaaatggaatccattgggtggtggtatggccgcacagggggcatggtggtcatcacgggtgctgctctgtgcctatgttctgctggtgctgggtctcatctgcccctgacactgggtcc
cgcaagtaa