

IL17B-interleukin 17B Gene

Reference: PTXBC113925

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IL-17B; IL-20; NIRF; ZCYTO7; interleukin-17B; cytokine-like protein ZCYTO7; interleukin 20; interleukin-17 beta; neuronal interleukin-17 related factor; interleukin 17B		

For research use only.

GENE DESCRIPTION	interleukin 17B
GENE ID	27190
ORIGIN SPECIES	Human
ACCESSIONS	BC113925
VECTOR	pDONR223
NCBI SYMBOL	IL17B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactggcctcacaacctgctgtttcttaccattccatctctctggggctgggcccagggagcccaaaagcaagaggaggggcgaaggcgctgggcccctggcccctggcccctcaccaggtgccactggactggtgtcacggatgaaacctatgccgcatgga
ggagtatgagaggaaacatcgaggagatggtggcccagctgaggaaacagctcagagctggcccagagaaagtgtgaggtcaacttgcagctgtggatgtccaacaaggagagcctgtctccctggggctacagcatcaaccagacccagccgtatccccgtgacctgccggagg
cacggctgctgtgtctgggctgtgtgaaccccttcacatgcaggaggaccgcagcatggtgagcgtgccggtgttcagccaggttctctgtgcgccgcccctctgcccgccaccgccgcacaggccttggccgacgcgcagtcattgagaccatcgtgtgggctgcacctg
catctctga