



Reference: PTXBC114429

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIF3; EIF3S10; P167; TIF32; eIF3-p170; eIF3-theta; p180; p185; eukaryotic translation initiation factor 3 subunit A; EIF3, p180 subunit; centrosomin homolog; cytoplasmic protein p167; eIF-3-theta; eIF3 p167; eIF3 p180; eIF3 p185; eukaryotic translation in		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 3, subunit A
GENE ID	8661
ORIGIN SPECIES	Human
ACCESSIONS	BC114429
VECTOR	pDONR223
NCBI SYMBOL	EIF3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr