

## DOK1-docking protein 1, 62kDa (downstream of tyrosine kinase 1) Gene

Reference: PTXBC114440

|                               |                                                                  |                     |                |
|-------------------------------|------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                              | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | P62DOK; pp62; docking protein 1; Downstream of tyrosine kinase 1 |                     |                |
| <b>For research use only.</b> |                                                                  |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | docking protein 1, 62kDa (downstream of tyrosine kinase 1)                                       |
| <b>GENE ID</b>          | 1796                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC114440                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | DOK1                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggcagcagcagtgatggaaggccgttttttgcagagtcagcgtttgggaccaagagtgaggagagacgtggccgtgctctaccggccagtcgccagggcgtagcgggctcaggtctttaccataagggtgcagctctgggggtggccgaggagctcgcgcgc  
ctggactgcaaaagtgcctgctgctgagtggtgagtggtggcccgctaccgtggagacccccctgagccggcgccactgcttccgctggacactgctcagcgtcgcacgtgctggcgccgacgcgcgtccagtgacagctgggtgcagagctgtgcccgaacgcct  
ttccgaaggcagctgactctgctgctaccgataaccacctaagcttctgcccctggagatgctggagaactcctttgacagccctactgggaaggatcccaatctgggtaacggtgcagaggactgagggcccgagcgtgtggcctgcatggtcctactgtgaggtgg  
aggctgaaggctgactctcctgaccgtggggggccagagtcagatactggagccactcctgctcgtggccctacactctgttgcgtcgtatggccgggacaaggctcatgttctcttgcaggccggcgccgctgcccctcaggccctggaaccttcacettccagagggcacaggaa  
atgacatctccaggcagtgagactgcccaccggcagaaggccagggaaggccggcagaggggcagatgttctcagagctgactccatgaaggggaggtggcagaggggaagtggccttcccacctggccccaagagctcctgcagctccccagccctgtatgctg  
agcccttagactccctgctgctgctccatgcccctccaggagctcctatactacagaccccttgacagcacgtctgctcaggcaggagaggagtagacaacggaagaaacctctctattgggactgtatgagcatgctgcagcagcaggttgcgaaggccaagctgacagacccaag  
aggatcccatctatgaaacctgagggcctggcccgctccctccagggcctttatgactgctcgggagcccaaggatgcatggtggtgccaagctcgggtgaaaggaggaggggtatgagctccctacaacctgcccactgatgactacgtgtgcccacctcggagcaca  
aagccccctctgctcccaagccccaggcccgctccctgaacctggtactgcaactggcagtggtgcatcaaaagccacaactcagccctgtacagcaggttcagaagagcggggcctcaggagctgggactgtggctctctagtagtaggactgacaagactgggtcaagt  
cagaggctctactga