

CRYGA-crystallin, gamma A Gene

Reference: PTXBC114456

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRY-g-A; CRYG1; CRYG5; gamma-crystallin A; crystallin, gamma 1; gamma-crystallin 5; crystallin gamma A		

For research use only.

GENE DESCRIPTION	crystallin, gamma A
GENE ID	1418
ORIGIN SPECIES	Human
ACCESSIONS	BC114456
VECTOR	pDONR223
NCBI SYMBOL	CRYGA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcgtcccaattaccaggccaccagtacttctgcgccgagaccagctgcacaagttaaggctgtacgagagagatgactaccgagccttatgtctgagctcactgatgactgcgcctgtgtccagaactgttcctccctgagatctattccctccacgtactggaggctgct
gggtcctctatgaaatgcccaactaccggggcgcgagtatctgtgagcctgggactacagaaggtaccacagctggggggtgcagatgccaaagtcggctctttgagacgggtcaccgattttgactaa