

## SLAMF6-SLAM family member 6 Gene

Reference: PTXBC114495

|                |                                                                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CD352; KALI; KALIB; Ly108; NTB-A; NTBA; SF2000; SLAM family member 6; NK-T-B-antigen; NTBA receptor; activating NK receptor; natural killer-, T- and B-cell antigen |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | SLAM family member 6                                                                              |
| <b>GENE ID</b>          | 114836                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC114495                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | SLAMF6                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgttggtgctgttccaatcgtcctgtttgtctctgtttggcccagggaatgtatgttcacaaagcagcttaacccattgatgggaacgggattctggggagtcagtaactctccctggagtttctgcaggagagaaggtcaacttcacttggctttcaatgaaacatcttgccttcatagtaccccatgaaacaaaagtcagaaatccacgtgactaatccgaacagggaaagcgactgaacttcacccagtcctactcctgcaactcagcaactgaagatggaagacacaggtctttacagagccagatatccacaaagacgtctgcaagctgtccagttacacttgaggatattagacaactgaggaacatacaagttacaaatcacagtcagctatttcagaatatgacctgtgagctccatctgacttgcctctgtggaggatgcagatgacaatgtctcattcagatgggagggccttgggaacacacttcaagtcagccaaacctcactgtctctgggaccccaaggattccagtgaacagactacacctgcatagcagagaatgctgtcagtaattatccttctgtctctgccagaagctttgcgaagatgttaaaattcaatatacagataccaaatgattctgtttatgtttctgggatatgcatagtctcgtttcatcactgtctgttactgttttgaggaaagaagagattccctatctttgtctactcagcgaacacagggcccccagagatccgcaaggaacctagagatgtttcagtgctcccaacgaacacactgtgtatgcttcagtcactcattcaaacagggaacagaaatctggacacctagagaaatgatactatcacaatttactccacaattaatcattcaagagagtaaacccattttccagggaactgcccttgacaatgtcgtgtaa