



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C1orf215-chromosome 1 open reading frame 215 Gene

Reference: PTXBC114502

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C1orf215; C1orf215; LAP18; Lag; OP18; PP17; PP19; PR22; SMN; stathmin; leukemia-associated phosphoprotein p18; metablastin; oncoprotein 18; phosphoprotein 19; phosphoprotein p19; prosolin; stathmin 1/oncoprotein 18; testicular tissue		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 215
GENE ID	149421
ORIGIN SPECIES	Human
ACCESSIONS	BC114502
VECTOR	pDONR223
NCBI SYMBOL	C1orf215
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr