

FTSJ2-FtsJ homolog 2 (E. coli) Gene

Reference: PTXBC114514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FTSJ2; FJH1; HEL97; RRMJ2; rRNA methyltransferase 2, mitochondrial; 16S rRNA (uridine(1369)-2'-O)-methyltransferase; 16S rRNA [Um1369] 2'-O-methyltransferase; FtsJ RNA methyltransferase homolog 2; FtsJ homolog 2; MRM2 RNA methyltransferase homolog; cell d		

For research use only.

GENE DESCRIPTION	FtsJ homolog 2 (E. coli)
GENE ID	29960
ORIGIN SPECIES	Human
ACCESSIONS	BC114514
VECTOR	pDONR223
NCBI SYMBOL	FTSJ2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtacttgaagctggtgtgtttcttcagcgtcaagggttccacactgttgggagtcgctgcaagaatcggacagggcgtgagcactgtggctgacccgacatctcaggacccattgtgaaggctgcgaaggtggagagttaccgggtgcgaagcgcctcaagctcct
ggaggtgaacgagagccaccagattctcgcgcccagccttcgggtgttagactgtggggcagctcctggggcctggagtcaggtggcgggtgcagaaggtcaacgccgagcagacatccagctcctctgttggcttcgtgcttgggttagatcttcttcacatattccccctggaagg
agcaactttctgtccctgctgactgactgacccgagaacctcagagaaatctcagagtgcttctggcaggagagcagatgtattctgagcgacatggcgcccaatgccacaggggttcgggacctcgtatcagaggtcctcagcctgtgctgaccccttcagcgtgacc
ccagacatctgcaacctgggggacattccttgaataacctgggctgggaagtcgaagccgtcgggttacagaggagactgacagaggaattcagaatgtaaggatcatcaaacctgaagccagcaggaagagtcagaaagtgtacttcttggccacacagaccacgaaggaa
gggcactgtgaagcagtga