

KLRB1-killer cell lectin-like receptor subfamily B, member 1 Gene

Reference: PTXBC114516

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD161; CLEC5B; NKR; NKR-P1; NKR-P1A; NKR-P1A; hNKR-P1A; killer cell lectin-like receptor subfamily B member 1; C-type lectin domain family 5 member B; killer cell lectin-like receptor subfamily B, member 1; natural killer cell surface protein P1A; killer c		

For research use only.

GENE DESCRIPTION	killer cell lectin-like receptor subfamily B, member 1
GENE ID	3820
ORIGIN SPECIES	Human
ACCESSIONS	BC114516
VECTOR	pDONR223
NCBI SYMBOL	KLRB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaccaacaagcaatatatgctgagttaaactaccacagactcaggccagaaagtcttcacctcatctctcctcgggatgtctgtcaggggttcaccttgccatcaattgccctgaaactagctgtgctgggattattctcctgtcttgggttactgggtgagtggttcagtgacatc
cttaatacagaaatcatcaatagaaaaatgcagtggtgacattcaacagagcagggaataaaacaacagagagaccgggtctcttaaaactgcccaataatattggcagcaactccgagagaaatgctgttattttctcacactgtcaacccttggataaacagctagctgattgtccaccaaag
aatccagcctgctgttattcgagataaggatgaattgatacacacagaaacctgatacgtgacaagcaattctgttttgattggattaaattttcattatcagaaaagactggaagtggataaacggctctttttaattctaatgacttagaaattagaggtgatgctaagaaaacagct
gtatttccatctcacagacatctgttattctgagtactgtatagacagaaatcagatggatctgccaaaagaactaacctgtgagaataaagtgtatcctgactcttga