

MYH1-myosin, heavy chain 1, skeletal muscle, adult Gene

Reference: PTXBC114545

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL71; MYHSA1; MYHa; MyHC-2X/D; MyHC-2x; myosin-1; epididymis luminal protein 71; myHC-IIx/d; myosin heavy chain 2x; myosin heavy chain IIx/d; myosin heavy chain, skeletal muscle, adult 1; myosin, heavy chain 1, skeletal muscle, adult; myosin, heavy polyp		

For research use only.

GENE DESCRIPTION	myosin, heavy chain 1, skeletal muscle, adult
GENE ID	4619
ORIGIN SPECIES	Human
ACCESSIONS	BC114545
VECTOR	pDONR223
NCBI SYMBOL	MYH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr