



Reference: PTXBC114547

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cylicin-2; cylicin II; cylicin, basic protein of sperm head cytoskeleton 2; multiple-band polypeptide II; cylicin 2		

For research use only.

GENE DESCRIPTION	cylicin, basic protein of sperm head cytoskeleton 2
GENE ID	1539
ORIGIN SPECIES	Human
ACCESSIONS	BC114547
VECTOR	pDONR223
NCBI SYMBOL	CYLC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1