

C18orf45-chromosome 18 open reading frame 45 Gene

Reference: PTXBC114561

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms:	transmembrane protein C18orf45; C18orf45; hVVT; transmembrane protein 241		

For research use only.

GENE DESCRIPTION	chromosome 18 open reading frame 45
GENE ID	85019
ORIGIN SPECIES	Human
ACCESSIONS	BC114561
VECTOR	pDONR223
NCBI SYMBOL	C18orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgtcgtgaggagctctctggtcggccacctttgtacctgctacctggtctctfacctcagcaacaagtatgtgctgtctgttgaattfacctacctacattattccaagggtggcagacgctcattggtgactttgcttcatgtctctgaaactgggctggtagatcaacagc
agttcaagatctatgttctgtggtctctctcagtcctgtttgtgggtataatctatgctgggtccagagcattgtccagactggccattcctgtgtttctactttgcataatgtagctgaagttatcatctgtgggtaccagaagtgtttcagaagagaaaacatctcctgcaaatctct
agtcctctctctctctgcccagcagcagcagccttccctcaatgactcccagtttaaccagatggatattttggctataatcacttactctgtgtagggcctataaaaattctacagaagtcccagaaccagtgcttaagtgacattgaccagcaatactaaactatattcagttgg
tgctctggcatttgcactcaccacaggtgafctctcagcgtcctgacttcccattctgtacttctacagattccatggtagctgtgccaaggattttgggatttctctatgtttcagttacagtgaaagctaaaaaccctctggccccaggcagtgtagcagcctgatttcttgcct
agataatcacagctgcttatcaatattgctgttgatgcatcctgaccagtgcaaccacgggatgctcctctgctcgttgcgcttgagaggccttctgttttctcagagcgaagagctcctga
```