

PPL-periplakin Gene

Reference: PTXBC114620

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	190 kDa paraneoplastic pemphigus antigen; 195 kDa cornified envelope precursor protein		

For research use only.

GENE DESCRIPTION	periplakin
GENE ID	5493
ORIGIN SPECIES	Human
ACCESSIONS	BC114620
VECTOR	pDONR223
NCBI SYMBOL	PPL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr