

SNRNP48-small nuclear ribonucleoprotein 48kDa (U11/U12) Gene

Reference: PTXBC114630

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf151; dJ336K20B.1; dJ512B11.2; U11/U12 small nuclear ribonucleoprotein 48 kDa protein; U11/U12 snRNP 48 kDa protein; U11/U12 snRNP 48K; U11/U12-48K; small nuclear ribonucleoprotein 48kDa (U11/U12); small nuclear ribonucleoprotein, U11/U12 48KDa subunit		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein 48kDa (U11/U12)
GENE ID	154007
ORIGIN SPECIES	Human
ACCESSIONS	BC114630
VECTOR	pDONR223
NCBI SYMBOL	SNRNP48
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggcgagcctccacctgtggaggagcggcggcgctgcaggagagctgaacgagttcgtggagagcggcgtccgacgttgaggaggtgaccgcgtccctggcgtgggacctagatagcttgatctcggggaagaggaggcggcgaggatgaagttgtgatgt
ccatagcattccaatcatcacatgcctaaatcatcttggcaagcacatggcatctgtgattgaggaaaatgggtataccaagaagaagagatgaaatgtataatcctgagttttctatgaaaatgtgaagataccttgatttgaataaggactcacaatccagataattaaacaa
gctagaactgcagttgggaagacagtgatttataatcaagaattattcttcattgctgtggaagttcctttgaatcacaaacgggtttgttgatcactcaagctgacgtcttgcctctatgatttctgattgaggagacaaagaaagcgcctctgattctcaaatatgaaatgac
agcgtatctttagacttgctgccaaatcaatcaagacaatagtcgaaaaagtcgaaatcctgaaatcctggcagaagtcagagattataaaagagacgccagtcctatagagccaagaatgttcacataaccaagaatcatactgaggtgattcgagatgtgataatgt
gcacatggaagaactcagcaatcattggcagaagagcagagaagcagagggatgatgccgaaaagaatgaagaagcgcagctcagtcagttcagtcagtcgtggaagctatttggatgctgagtttcacgacatagaaggataggagtagaagccacataaaga
aaaagaacaaagataaggataaaaactgtgagtcgagaagaaggaaagaggatggggaagacacacatagtcataaaagaagaagcaaaaataaa