



Reference: PTXBC114924

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SAHH; adoHcyase; adenosylhomocysteinase; S-adenosyl-L-homocysteine hydrolase; S-adenosylhomocysteine hydrolase		

For research use only.

GENE DESCRIPTION	S-adenosylhomocysteine hydrolase
GENE ID	191
ORIGIN SPECIES	Human
ACCESSIONS	BC114924
VECTOR	pDONR223
NCBI SYMBOL	AHCY
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1