

## SLC35F1-solute carrier family 35, member F1 Gene

Reference: PTXBC114925

|                               |                                                         |                     |                |
|-------------------------------|---------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                     | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                    | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | C6orf169; dJ230I3.1; solute carrier family 35 member F1 |                     |                |
| <b>For research use only.</b> |                                                         |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | solute carrier family 35, member F1                                                              |
| <b>GENE ID</b>          | 222553                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC114925                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | SLC35F1                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgatccccctgagcagccgcagcagcagctgcagccgctgcgagcccccgcgaacatgtgtgaccaccatcgagaacctgcggccgaggcgagcgccgagccggcgagcctgtccgctctcccggtggcggtgcgagagggatccgcaagtgcgaacag  
ggagatgttaactctgtgcccctagccaggtgtatccctcttattgtggaattggctgactagcaagtatctgtcagaagattccagcccaacacacagctctccagagttctcaattacttctcttctgtgtctataccacacactagccgtcagacaaggagaagaaacct  
cctggcaattttacgacgaagatggtggaagtacatgattttgggactcatagacctggaagcaattatctgtgtgctcaaggcttaccatacacaactctgaccagtatccagctcctggactgttttggatcccgatcgtgattttgctctcctgtgtctcctgctgacccggtacaagcgtgt  
gcatttcatcggcatcgtgtgtcgcacatcctgggaatgggctgcatgtgtggagcagatgtgcttgtgggaagacatcagggagcagggggaataaagctggttaggggacacttctgtgtcttagggaggagccacactctatggtatttcaacgtctctgggaagaatcatcaccgaactctgag  
ccgagtggaattcctgggaatgattggtctcttggagcattttcagtggaattcaattagctataatggagcataaaggaaactgttgaaggtgcccctgggactggcaaataggactgctctacgttggcttttagtgcctgcatgtttgtgtctctacagctttatgccagtcgtcataaagaaaacca  
gtgccactcagtcacactcctctgctcacagcagactgtacagcctgttctgtgattgtttctctccactacaagtttcaggactttatctcctgtcttcttccacatcctcattgggctgtgtctctactcctccacctccacatagcccaggacccccgagtgataagcagttccgca  
atccttcaggacctgtgttggaattaccgaccacagctcaggtggaacctcagtcacctacaccagcctgggcccaggagaccgaagaggagcctcatgttctgtgtgcttag