

## C20orf152-chromosome 20 open reading frame 152 Gene

Reference: PTXBC114945

|                |                                                                                                                                                                               |                     |                |
|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                           | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                          | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C20orf152; CNMPD1; cyclic nucleotide-binding domain-containing protein 2; cyclic nucleotide (cNMP) binding domain containing 1; cyclic nucleotide binding domain containing 2 |                     |                |

**For research use only.**

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 20 open reading frame 152                                                              |
| <b>GENE ID</b>          | 140894                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC114945                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C20orf152                                                                                         |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaggagacatatgtaacttatgcctggcagctcctgaagaaggaaactgggactgtaccagctgccatggatatcatcataatgatccgagtggtgtaaaatgttcgccaaggcctcaggggattccgggaatatcaaatcattgagactgctcactggaagcacctatctctctctg  
ggataaaaagatgcaaacgagtcacattgataccatggactcattgcagaggagggtcactttctccaaggccattcagatcatgcagaagaagccttcctggagaacagagatgagatccaggccgtctgttaacatcttcaggtttctggatagctatcggaactacgcagagc  
ccctgcagctgctctggccaaagtcacgtttgaacggtttgtgcagggcgtgtgatcatcaagaaggggcagaaggggcaacagcttttattcatctacctgggcacagttgcaataaccaaggacgagatggcagcagtgctctctagatccccaccgaaattgctgcacaag  
ggtagctgttttgggaaatggacgtttctgcatgcttcagtgaggaggtccaccatcgtctgtatggaagaacggagttcctggtttgtgacccgggaggactctttgctaataagctggaccaggaagttcagaagatgctcagatcgtttgaaatttttaggaagatggagctgtttgcat  
catggtctgatgagaagctctggcagctggttagccatggcgaagatagagaggttctgtatggcagctgactcaaaagatttggagagtcaccctcatcatgtttatcagcaaggcgagctgtgaagtcctgctggcgtgttgacacttggggcctcccttccaccgtagatggatcgt  
gcagcacctggagctgatagatggcagacctctgaagaccacctgagtgaaactctctatggaaagatttaagggaattccagatcaaatcatatcctctgcaagacttttagctccttgaactccacatctcaaaaaagcctgggggctacaggggacaagcttcagcaggaagatcag  
aacctcaggagacactctcccaagatgctggggccgaagatccaatccaggcctcgtcagtcgatcaaatgtgccatgatcaatatcaagcctggtgagctcccaaggagctgcagctgggggctacgtgaaggtgcacactgtggagcagggagaaattttgtga