

ASB14-ankyrin repeat and SOCS box-containing 14 Gene

Reference: PTXBC114953

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ankyrin repeat and SOCS box protein 14; ankyrin repeat domain-containing SOCS box protein Asb-14; ankyrin repeat and SOCS box containing 14		

For research use only.

GENE DESCRIPTION	ankyrin repeat and SOCS box-containing 14
GENE ID	142686
ORIGIN SPECIES	Human
ACCESSIONS	BC114953
VECTOR	pDONR223
NCBI SYMBOL	ASB14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgaaatgcatttgaattctttcagagctctaaagatactgattccagttacggatcttgccttaagcagagtgaggatcagtcagttcactgtcagcagcaggagcacatcctcagtcgctggaactcctcatccaggctgatttgatgtgaactcatgctggatcagagaat
taacaaactacgatgaccacagggaagtcagcttctgattttgctgtatcaaacagtgacctctctcagtaagctgcttctgagtgctggagctctgcctaatcaagaccgggttaactgcctccagatagccctcaggatgggcaactatgagctgacagctgctgctaaggcatggg
ccaatgtcaattactctgcagagtttaacctttacatttccatcagcactgcaatacactctgaaagatgaagtcagctcaggatgctgctgaactatgggtatgacacagagcagatgttttgattgcccacatggagacaagtcacatcctctatactgttgaaggctggacatctacagtt
atcaagatactaaagtaagtcacagattggctgtattcccttacagtcatttttggccctgggaactatggcaataatagattacatccactctcagaggcttggttgaagagatcctag