



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

DTWD2-DTW domain containing 2 Gene

Reference: PTXBC114955

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DTW domain-containing protein 2; DTW domain containing 2		
For research use only.			

GENE DESCRIPTION	DTW domain containing 2
GENE ID	285605
ORIGIN SPECIES	Human
ACCESSIONS	BC114955
VECTOR	pDONR223
NCBI SYMBOL	DTWD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagatcgcagaagaaggcacgaacactccaggagccggtgtgcgcgctcttggggccctcaagctctcagacgcgaacacagaggagcggcgaggagggcgcgagfgcggcgccgctgctctggcgacagaggcgacgacagatgcggacgggctgtggag
ctgcgggtggagcggcgagcggaggcctgagtgccaccgctgcagcggcctcagaagtgtgtgtgtcctattctccagcgcacctctgcatatctcaccactgtacataattcagcatccagcagaggaaaaaagaagtgtgtcgtacgtctctactacgacgatgcctc
cccagacagatgtaaaagtgaagatcggctgtcttcagtgaaagaagaatcattctcaactgtttggcagcttgcgtacattatgatattatccagggcgctgaagctcctaatttgaagaatttatattagattctctgttacttaccatcatcatcattgatgtacatgc
gcagggctaaaggctcttctaagaagctcgtgtcttcgacatcccaaaacagctgtgccaatttcaagcttctagtcagatgtgaattcgaatgcagcttctcactgcggatgtgcgctgtgtcctcttcattctggcataaaataattacacagaagacttgcct
ctgcctctcaagctttatgctcttcttcacagatcagatgtgcgccaattgcctcgaagaagaacctctgaaaattgattatctcttaaccacaaatgcgaagaagaacaaacgaacactcgaagaatgcgaatgttaataagatgtgttaaaattaaat