

## RSPO1-R-spondin homolog (*Xenopus laevis*) Gene

Reference: PTXBC114966

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CRISTIN3; RSPO; R-spondin-1; R-spondin homolog; roof plate-specific spondin-1; R-spondin 1		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	R-spondin homolog ( <i>Xenopus laevis</i> )
<b>GENE ID</b>	284654
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC114966
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	RSPO1
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggcctgggctgtgtgtggccctgttctgagctggacgcacctcaccatcagcagccgggggatcaagggaaaaggcagagcggatcagtcggagggagccagcctgtgccaaaggctgtgagctctgctctgaagtcaacggctgctcaagtgtctaccaagctgttcacctgtgagaggaacacatccgaggtggcgtctgtcctgcccactggatactcgacgccgcaaccccgacatgaacaagtcacataaatgcaagatcgagcactgtgagcctgttcagccataactctgcacaaagttaaggaggctgtacctgcacaagggcctgtctatccagctgtcccgaggctctcagctgccaatggcaccatggagtgagtagtctcgcgcaatgtaaatgagcagtggtctcctggggccctgtccaagaagcagcagctctgtgttccggagggtccgaggagcggacacgcaggggtgctacatgccctgtggggaccatgctgctctgacaccaaggagaccggaggtgcacagtgaggagagtgccgtgtctgaggggcaagagaggaggaaggagggccaggccggcgggagaatgccaacaggaacctggccaggaaggagcaaggaggggggtgctgtctcgaagacgcaagggcagcaacagcagcagcagcaagggacagtggggccactcacatctcagggcctgctag