



Reference: PTXBC115039

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BSAC; MAL; MRTF-A; MKL/myocardin-like protein 1; basic, SAP and coiled-coil domain; megakaryoblastic leukemia 1 protein; megakaryocytic acute leukemia protein; myocardin-related transcription factor A; megakaryoblastic leukemia (translocation) 1		

For research use only.

GENE DESCRIPTION	megakaryoblastic leukemia (translocation) 1
GENE ID	57591
ORIGIN SPECIES	Human
ACCESSIONS	BC115039
VECTOR	pDONR223
NCBI SYMBOL	MKL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]