

TAS2R46-taste receptor, type 2, member 46 Gene

Reference: PTXBC115401

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	T2R46; T2R54; taste receptor type 2 member 46; taste receptor type 2 member 54; taste receptor, type 2, member 46; taste 2 receptor member 46		

For research use only.

GENE DESCRIPTION	taste receptor, type 2, member 46
GENE ID	259292
ORIGIN SPECIES	Human
ACCESSIONS	BC115401
VECTOR	pDONR223
NCBI SYMBOL	TAS2R46
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataacttttctgccatctttttccattctaatagtggtacatttggttggaattttgctaatggctcatagcattgtaattccattgagtggttcaagacaaaagatcttttctgaccaaattctcactgctctggcagctccagagttggtttactctgggtattatataaattg
gtatgcaactgagttgaatccagcttttaacagtatagaagtaagaattactgcttaaatgctgggcagtaatacaaccttcagcaactggctgctactagcctcagcatatttttctcaagattgccaatttcaccaactattttctcacttaaaaggagagtaagagttgtgtctg
gtgatactattgggctttgctatttttgggttgcatttttggataaacatgaatcagattatggacaaaagaataatgaaggaaacatgacttggaaatcaaacagagagtgcaatgtacctttcaaatcaaacggtaaccatcctagcaaaccttagttcccttcactctgacctgataatct
tttctgctttaactctgtctgtgtaaacatctcaaaaagatgcagctccatggcaaggatctcaagatccagcatgaaggtccacataaaagctatgcaaacctgtgacctcttctctgttatgtgccatttacttctgtccataatcatgtcagtttag