

MYOM1-myomesin 1, 185kDa Gene

Reference: PTXBC116183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SKELEMIN; myomesin-1; 190 kDa connectin-associated protein; 190 kDa titin-associated protein; EH-myomesin; myomesin (M-protein) 1 (190kD); myomesin 1 (skelemin) 185kDa; myomesin 1 variant 1; myomesin 1, 185kDa; myomesin family member 1; myomesin 1		

For research use only.

GENE DESCRIPTION	myomesin 1, 185kDa
GENE ID	8736
ORIGIN SPECIES	Human
ACCESSIONS	BC116183
VECTOR	pDONR223
NCBI SYMBOL	MYOM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr