

C1orf31-chromosome 1 open reading frame 31 Gene

Reference: PTXBC116455

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf31; CEMCOX4; cytochrome c oxidase assembly factor 6 homolog; cytochrome c oxidase assembly factor 6		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 31
GENE ID	388753
ORIGIN SPECIES	Human
ACCESSIONS	BC116455
VECTOR	pDONR223
NCBI SYMBOL	C1orf31
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggccggggaggtccctactgtcccgagccgggggttcctctgtgcaaaacgggggtggcactccaatgcctgcttgattgtggcccccacacctgtttctacagcgcttagctcatcgagtaggaatggcagccccctatgaaggaaagacaggtctgctggggg
ccgggatgagtactggaagtgtttagatgagaacttagagatgcttctcaatgcaagaagttaagaagctctttcgaatcaagtgtgccccacagtgataaaatattttgataaaagaagagactttaaaatcaagaaaaatttgaagcaggacaatttgagccttcagaacaactg
caaatcctag