



Reference: PTXBC117116

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PAC1; PAC1R; PACAPR; PACAPRI; pituitary adenylate cyclase-activating polypeptide type I receptor; PACAP receptor 1; PACAP type I receptor; PACAP-R1; adenylate cyclase activating polypeptide 1 (pituitary) receptor type I; pituitary adenylate cyclase activa		

For research use only.

GENE DESCRIPTION	adenylate cyclase activating polypeptide 1 (pituitary) receptor type I
GENE ID	117
ORIGIN SPECIES	Human
ACCESSIONS	BC117116
VECTOR	pDONR223
NCBI SYMBOL	ADCYAP1R1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]