



Reference: PTXBC117117

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMIT; proton myo-inositol cotransporter; H(+)-myo-inositol symporter; h(+)-myo-inositol cotransporter; proton (H+) myo-inositol symporter; solute carrier family 2 (facilitated glucose transporter), member 13; solute carrier family 2 member 13		

For research use only.

GENE DESCRIPTION	solute carrier family 2 (facilitated glucose transporter), member 13
GENE ID	114134
ORIGIN SPECIES	Human
ACCESSIONS	BC117117
VECTOR	pDONR223
NCBI SYMBOL	SLC2A13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr