



Reference: PTXBC117131

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPND; FESD; GluN2A; LKS; NMDAR2A; NR2A; glutamate receptor ionotropic, NMDA 2A; N-methyl D-aspartate receptor subtype 2A; N-methyl-D-aspartate receptor channel, subunit epsilon-1; N-methyl-D-aspartate receptor subunit 2A; glutamate receptor, ionotropic, N		

For research use only.

GENE DESCRIPTION	glutamate receptor, ionotropic, N-methyl D-aspartate 2A
GENE ID	2903
ORIGIN SPECIES	Human
ACCESSIONS	BC117131
VECTOR	pDONR223
NCBI SYMBOL	GRIN2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1