

KTN1-kinectin 1 (kinesin receptor) Gene

Reference: PTXBC117132

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CG1; KNT; MU-RMS-40.19; kinectin; CG-1 antigen; kinesin receptor; kinectin 1		
For research use only.			

GENE DESCRIPTION	kinectin 1 (kinesin receptor)
GENE ID	3895
ORIGIN SPECIES	Human
ACCESSIONS	BC117132
VECTOR	pDONR223
NCBI SYMBOL	KTN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr