



Reference: PTXBC117137

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POLS; LAK-1; LAK1; POLK; TRF4-1; TRF41; TUTASE5; non-canonical poly(A) RNA polymerase PAPD7; DNA polymerase kappa; DNA polymerase sigma; PAP associated domain containing 7; PAP-associated domain-containing protein 7; TRAMP-like complex polyadenylate polym		

GENE DESCRIPTION	polymerase (DNA directed) sigma
GENE ID	11044
ORIGIN SPECIES	Human
ACCESSIONS	BC117137
VECTOR	pDONR223
NCBI SYMBOL	POLS
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]