



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LY6K-lymphocyte antigen 6 complex, locus K Gene

Reference: PTXBC117142

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CT97; HSJ001348; URLC10; ly-6K; lymphocyte antigen 6K; cancer/testis antigen 97; up-regulated in lung cancer 10; lymphocyte antigen 6 complex, locus K

For research use only.

GENE DESCRIPTION	lymphocyte antigen 6 complex, locus K
GENE ID	54742
ORIGIN SPECIES	Human
ACCESSIONS	BC117142
VECTOR	pDONR223
NCBI SYMBOL	LY6K
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]