



Reference: PTXBC117145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EHB3; PRP17; PRPF17; pre-mRNA-processing factor 17; EH-binding protein 3; PRP17 homolog; cell division cycle 40 homolog; hPRP17; pre-mRNA splicing factor 17; cell division cycle 40		

For research use only.

GENE DESCRIPTION	cell division cycle 40 homolog (S. cerevisiae)
GENE ID	51362
ORIGIN SPECIES	Human
ACCESSIONS	BC117145
VECTOR	pDONR223
NCBI SYMBOL	CDC40
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG