

ST18-suppression of tumorigenicity 18 (breast carcinoma) (zinc finger protein) Gene

Reference: PTXBC117147

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ST18, C2H2C-type zinc finger; NZF3; ZC2H2C3; ZC2HC10; ZNF387; suppression of tumorigenicity 18 protein; neural zinc finger transcription factor 3; suppression of tumorigenicity 18 (breast carcinoma) (zinc finger protein); suppression of tumorigenicity 18,		

For research use only.

GENE DESCRIPTION	suppression of tumorigenicity 18 (breast carcinoma) (zinc finger protein)
GENE ID	9705
ORIGIN SPECIES	Human
ACCESSIONS	BC117147
VECTOR	pDONR223
NCBI SYMBOL	ST18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1